


## ACCESSORIES



### INFUSION PUMP HOLDER

Holds up to 3 pumps  
P-1080-G Infusion Pump Holder



### CV HOOK

Molded polymer six hook IV assembly  
P-1090-CV-TH Assembly



### 8-HOOK ASSEMBLY

Extended IV Hook Assembly  
P-1080-H Stainless Steel



### S'MORE HOOK

For mounting additional IV bags  
P-1080-J Stainless Steel



### FOLEY HOOK

P-1080-D Foley 2 Hooks  
P-1080-E Foley 1 Hook  
(Factory installed)



### IV SET ORGANIZER

Also holds secondary bags  
P-1080-B Stainless Steel



### PATIENT HANDLE

Adjustable handle provides easy mobility.  
P-3301 Stainless Steel  
P-3302-BLK Black  
P-3303 CuVerro®



### WIRE BASKET

P-3201 Large (Stainless)  
P-3202 Large (Black)  
P-3211 Med (Stainless)  
P-3212 Med Black



### O<sub>2</sub> HOLDER

Holds E-size oxygen tanks  
P-3601 Stainless Steel



### OUTLET STRIP

Hospital grade  
P-3401 Six outlet strip



### INFUSION PUMP HOOP

Holds up to 6 pumps  
P-1080-F Stainless Steel  
P-1083-F CuVerro®



### MULTIPLE PUMP HOLDER

Infusion pump holder  
P-1080-C Stainless Steel  
P-1083-C CuVerro®



infusion pump stands & IV poles



Copper is an intrinsic bactericidal element and **CuVerro®** is the **ONLY** solid touch surface metal registered by the EPA to continuously kill bacteria<sup>1</sup> that cause infections and pose a risk to human health. Lab tests show more than 99.9% of bacteria<sup>1</sup> associated with disease and infection, including MRSA<sup>1</sup>, are killed within two hours when in contact with copper. Bacteria will continue to be killed 24 hours a day, week after week, between regular cleanings.

- MRSA
- Enterobacter aerogenes
- Pseudomonas aeruginosa
- Escherichia coli O157:H7
- VRE
- Staphylococcus aureus

**Pedigo Products EPA Establishment No. 088754-WA-001**

The use of CuVerro is a supplement to and not a substitute for standard infection control practices; users must continue to follow all current infection control practices, including those practices related to cleaning and disinfection of environmental surfaces. This surface has been shown to reduce microbial contamination, but it does not necessarily prevent cross contamination. CuVerro is a registered trademark of GBC Metals, LLC and is used with permission. PP-0004-1202





## infusion pump stands & IV poles



**P-1080**




**P-1083**



**P-1090-CV**



**P-1572**



**P-1576**

<b>Infusion Pumps</b>	Up to 6 pumps at 10 lbs each	Up to 6 pumps at 10 lbs each	Up to 4 pumps at 10 lbs each	Up to 2 pumps at 10 lbs each	Up to 2 pumps at 10 lbs each
<b>Total MAX Weight Capacity*</b>	125 lbs	125 lbs	90 lbs	70 lbs	70 lbs
<b>Base Weight</b>	27 lbs	27 lbs	19 lbs	11 lbs	11 lbs
<b>Base Material</b>	Powder coated steel	Powder coated steel	High impact polymer	Stainless steel	Stainless steel
<b>Footprint (outer perimeter)</b>	23"	23"	23"	21.25"	21.25"
<b>Interlocking Bases</b>	No	No	Yes (max of 6)	No	No
<b>Caster</b>	3" precision ball bearing - sealed	3" precision ball bearing - sealed	3 "Precision ball bearing - sealed	2" swivel dual wheel	2" swivel dual wheel
<b>Adjustable Height</b>	Hand Operated 64.75" - 124"	Hand Operated 64.75" - 124"	Hand Operated 74" - 122"	Hand Operated 67.5" - 117.5"	Foot Operated 67.5" - 117.5"
<b>Pole Material</b>	Stainless steel	CuVerro® bactericidal alloy	Stainless steel	Stainless steel	Stainless steel
<b>Bag Drop Dampening Mechanism</b>	Yes	Yes	Yes	No	No
<b>Hook options</b>	<ul style="list-style-type: none"> <li>• 6 hook stainless steel</li> <li>• 6 hook ClearView™ hook (CV)</li> <li>• 8 hook stainless steel (rake)</li> </ul>	<ul style="list-style-type: none"> <li>• 6 hook stainless steel</li> <li>• 6 hook ClearView™ hook (CV)</li> <li>• 8 hook stainless steel (rake)</li> </ul>	<ul style="list-style-type: none"> <li>• 6 hook stainless steel</li> <li>• 6 hook ClearView™ hook (CV)</li> <li>• 8 hook stainless steel (rake)</li> </ul>	<ul style="list-style-type: none"> <li>• 6 hook ClearView™ hook (CV)</li> <li>• 4 hook</li> </ul>	<ul style="list-style-type: none"> <li>• 6 hook ClearView hook™ (CV)</li> <li>• 4 hook</li> </ul>
<b>Accessories</b>	<ul style="list-style-type: none"> <li>• Patient Handle</li> <li>• Wire basket</li> <li>• O<sub>2</sub> holder</li> <li>• Outlet Strip</li> <li>• Foley hook</li> <li>• IV set Organizer</li> <li>• S'more Hook</li> <li>• Multiple Pump Holder</li> <li>• Infusion Pump Hoop</li> <li>• Infusion Pump Holder</li> </ul>	<ul style="list-style-type: none"> <li>• Patient Handle</li> <li>• Wire basket</li> <li>• O<sub>2</sub> holder</li> <li>• Outlet Strip</li> <li>• Foley hook</li> <li>• IV set Organizer</li> <li>• S'more Hook</li> <li>• Multiple Pump Holder</li> <li>• Infusion Pump Hoop</li> <li>• Infusion Pump Holder</li> </ul>	<ul style="list-style-type: none"> <li>• Patient Handle</li> <li>• Wire Basket</li> <li>• O<sub>2</sub> Holder</li> <li>• Outlet Strip</li> <li>• Foley Hook</li> <li>• IV set Organizer</li> <li>• S'more Hook</li> <li>• Infusion Pump Hoop</li> <li>• Infusion Pump Holder</li> </ul>	<ul style="list-style-type: none"> <li>• Patient Handle</li> <li>• Wire Basket</li> <li>• O<sub>2</sub> Holder</li> <li>• Outlet Strip</li> <li>• Foley Hook</li> <li>• IV set Organizer</li> <li>• S'more Hook</li> </ul>	<ul style="list-style-type: none"> <li>• Patient Handle</li> <li>• Wire Basket</li> <li>• O<sub>2</sub> Holder</li> <li>• Outlet Strip</li> <li>• Foley Hook</li> <li>• IV set Organizer</li> <li>• S'more Hook</li> </ul>

\* Applied at center of pole, 3 ft. MAX from ground.

