

Intelligent Fluid Warming

ivNext®



Entthermics is proud to offer a new IV Inventory Management System - the ivNext Fluid Warmer. Our WarmSafe™ Zone heating technology provides clinicians a simple way to monitor and precisely warm IV bags with a bag temperature tolerance of $\pm 0/-1.1^{\circ}\text{C}$ (-2°F).

- Compact size allows for placement on countertop or mobile cart for easy access to fluids.
- Warms 10 1-liter IV bags.
- Organized inventory allows for faster retrieval of fluids.
- Know exactly when an IV bag has reached the set point temperature - adjustable $32-40^{\circ}\text{C}$ ($90^{\circ}-104^{\circ}\text{F}$) injection mode.
- Visual indication of bag readiness.

Easy to Read

Easily recognizable symbols allow staff to act fast.



Blue = Warming



Green = Ready to Use



Yellow = Expired

Warming

Blue light

Visual indication when an iv bag has not yet reached set point temperature.

Ready to Use

Green light

Visual indication when an iv bag has reached set point temperature.

Expired

Yellow light

Visual indication when an iv bag has been warmed continuously for 14 days or longer.

