

Your ASC equipment solution!

✓ in stock ✓ fully assembled ✓ ready to use



CDS-3072-ADJ-E

Electric Adjustable Space Station

Space Station

Electric Adjustable Table

A specialty back table designed with additional room for large cases such as orthopedics, spinal fusions, neuro and craniotomies. The tuck-away cantilevered shelf holds multiple trays and is angled for clear observation of instruments.

Top Tier 150lbs (68kg)

Main Surface 200lbs (91kg)

Bottom Storage.. 200lbs (91kg)

LED Lights

- Easy to clean
- 6300K cool white light



new!



co-re caster™
collection



SG-80-A-SS
table



P-536-GS
seating



P-1120-L-SS
hamper

The **GUARDIAN**
with *CuVerro®* Side Rails



Guardian 7500 Series
transport / trauma stretcher



P-2022
blanket warmer

element™



P-1090-CV
iv stands

The use of CuVerro is a supplement to and not a substitute for standard infection control practices; users must continue to follow all current infection control practices, including those practices related to cleaning and disinfection of environmental surfaces. This surface has been shown to reduce microbial contamination, but it does not necessarily prevent cross contamination. CuVerro is a registered trademark of GBC Metals, LLC and is used with permission. PP-0004-1202